


©



Illustrating No. 2071 Trial Set, closed.
 Illustrating No. 2071 Trial Set, open.
 Illustrating No. 2080 Trial Set, open.

OFFICE STYLES, LEATHER COVERED—*Continued*


NUMBER	DESCRIPTION. 20 PAIRS SPHERICALS, 11 PAIRS CYLINDERS, CONTENTS "F." See Page 234
2070	Black seal grain leather covered, velvet lined case, removable tray. Lenses, 36 mm. in 1806 Nachet test rings.
2071	Same as No. 2070, except 1801 Alumnico test rings. Lenses, 36 mm.
2073	Same as No. 2070, except 1803 Alumnico test rings. Lenses, 32 mm.
2074	Same as No. 2070, except 1804 Nachet test rings. Lenses, 32 mm.
DESCRIPTION. 20 PAIRS SPHERICALS, 11 PAIRS CYLINDERS, CONTENTS "F." See Page 234	
2080	Black seal grain leather covered, velvet lined case, removable tray. Lenses, 36 mm. in 1806 Nachet test rings.
2081	Same as No. 2080, except 1801 Alumnico test rings. Lenses, 36 mm.
2083	Same as No. 2080, except 1803 Alumnico test rings. Lenses, 32 mm.
2084	Same as No. 2080, except 1804 Nachet test rings. Lenses, 32 mm.

No. 2080 style case is uniform in size with office cases designed to hold contents A, B, C, D and H
See illustration on opposite page.



Office Style Case Making

Ⓐ



TRIAL FRAMES

38 MM. NUMBER	33 MM. NUMBER
------------------	------------------

2387	2385	Patented double cell, Wells extension temples, celluloid scales on eyes and bar
2387 B	2383 B	Patented, same as No. 2387, except patent revolving cell and thumb screw
2387 C	2383 C	Patented, same as No. 2387, except patent revolving cell and thumb wheel
2467 B		Patented, same as No. 2387 B, except Rhoad's full circle scale and double axis pointer

- NOTE.—1. Above trial frames made with Wells patent temples unless otherwise ordered.
 2. Above trial frames made with three cells if desired, subject to an extra charge.
 3. All patent revolving cell trial frames may be had with lock as shown on frame 2397 B, page 269, if so desired, at a slight additional charge. Specify "with lock."
 4. Above trial frames can be made with D attachment for independent pupillary adjustment when so ordered, subject to a slight extra charge.
 5. Self-adjusting Bridge, shown on following page, is regularly supplied on above trial frames without extra charge.



TRIAL SETS, TRIAL FRAMES AND ACCESSORIES

NO branch of our business demands more careful treatment, more attention to detail than the manufacture of trial sets. The refractionist must largely depend upon the accuracy of his trial lenses for the correctness of his work, and in making the selection of this necessary adjunct to his equipment he has little to aid him, but must depend almost solely upon the standing and reliability of the maker.

It therefore behooves us to carefully guard the enviable reputation of our trial sets, which has already been instrumental in marketing our annual production of over 3000 sets.

The quality of stock, superiority of finish and durability of construction are details carefully considered in connection with accuracy of foci.

In connection with our comprehensive line of trial sets we would call especial attention to our various trial frames, many of the most desirable and practicable patterns of which are protected by Letters Patent owned by us. We believe we are justified in claiming this line includes the best and most popular trial frames made, which is best evidenced by the large quantity of these goods we are making.

INFORMATION FOR ORDERING

On the following pages will be found schedules showing the exact contents and list of foci of lenses regularly supplied in the various trial sets illustrated. We strongly urge the selection of some one of these various regular contents, as among them we feel sure an assortment will be found answering practically any requirement, and thereby delay in the filling of an order is largely avoided and expense saved.

All lenses, including prisms, are ground in the dioptral system, and focus strips are printed in this system only.

All trial sets except the pocket styles are regularly supplied with white celluloid strips, stamped with focus and filled with black, making them clear and indelible.

All of the larger styles of trial sets are made with removable trays, and, unless otherwise ordered, velvet lined trays are supplied. We would call especial attention to our all-wood trays, which are gaining in popularity owing to their cleanliness and durability.

All-wood Tray If an all-wood tray is ordered in connection with an all-wood case the order should be explicit. Thus: No. 2121 oak trial set with all-wood oak tray, or use No. 02121.

All trays, except in traveling style cases, are made skeleton construction, allowing dust to fall on removable pad under the tray, making it very convenient to keep the trial set clean.

All trays in office and dress suit style cases containing contents A, B, C, D, H are interchangeable.

Office Style All office style trial sets are covered with black seal grain leather and have velvet lined removable trays. If all-wood tray is desired order must so state. All cases have nickel-plated locks, catches, quadrant stops and hinges.

Office Style, Hard Wood Our oak cases are made of fully seasoned selected quartered stock, finely finished. All have removable trays velvet lined. If all-wood tray is desired order must so state. Cases are also furnished in mahogany if so ordered.

Automatic Tray Tilting Attachment This provides a means of automatically raising and lowering back of tray when case is opened or closed. This attachment will be furnished at a slight extra charge upon any office style cases except roll and glass front cabinets. See illustration, page 240.



Neutralizing Trial
Set Lenses



Tray Making

Office Style, Roll and Glass Front Cabinets Our upright cabinets are made of fully seasoned selected quartered oak, finely finished. All have removable trays velvet lined. If all-wood tray is desired order must so state. Trays are arranged to draw forward, making the lenses easily accessible. Cabinets can also be furnished in mahogany if so ordered.

Finish of Oak Cases Oak cases are finished antique color (light) unless otherwise ordered. Other finishes, as golden oak, mission, early English, Flemish, etc., will be furnished when order so specifies.



Traveling Style Case Making

Traveling, Dress Suit and Folding Style Cases are covered with black seal grain leather with nickel-plated trimmings, including quadrant stop for cover and with removable velvet lined tray with solid bottom. Space is provided underneath tray for stock, also with removable tray for interchangeable lenses, holding one gross pairs. The cover contains a pocket for holding test types. These cases are also furnished in tan leather if so ordered. Cases also furnished with oxidized bronze trimmings instead of nickel-plated if so ordered.

Sizes of Lenses and Trial Rings Trial lenses are made in two sizes and measure exactly 36.3 mm. and 31.5 mm. in diameter. For convenience they are designated as 36 mm. and 32 mm. respectively throughout this catalogue. These sizes are sometimes designated as XX (36.3 mm.) and XXX (31.5 mm.). See page 29, Introductory Section.

Trial rings are 37.6 mm. and 32.8 mm. outside and are designated as 38 mm. and 33 mm. for convenience. The metric system has been adopted with the publication of this edition, the inch system having been discontinued. The diameter of trial rings is denoted by their respective catalogue numbers. The unit figures 0 and 1 indicate the larger rings and unit figures 3 and 4 the smaller rings. It is, therefore, not necessary to state size of rings in ordering trial sets, the above being submitted merely as a matter of information.

We have excellent facilities for filling orders for trial sets on short notice, as a large number of the staple and some special styles are constantly kept in stock.

We are prepared to make any special style of trial set to order. It will be understood that orders for regular styles can be filled more promptly, as a trial set with special assortment of lenses or requiring a specially constructed case must of necessity be made from the beginning. Furthermore, the quantity of lenses or the assortment of foci, if differing even in a single pair from a regular assortment, cannot be changed without extra charge, as special focus strips or a specially constructed case and tray would be required, adding materially to the cost.